

# Purification Results from the LZ WIMP Search

**Carter Hall,  
University of Maryland**


**February 16, 2024**

**Nagoya Workshop on  
Technology and  
Instrumentation in  
Future Liquid Noble  
Gas Detectors**




# Experimental Strategy


- Online gettering of electronegatives in gas phase with SAES Megatorr / St707 technology
  - Continuous gas recirculation and purification at rates up to 500 SLPM driven by two all-metal compressors
  - Two-phase heat exchanger to mate the liquid and gas portions of the circuit
- Removal of trace of Krypton with gas charcoal chromatography
  - Performed off-site at SLAC prior to condensing into LZ cryostat
- Control of  $^{222}\text{Rn}$  through materials screening, cleaning, & etching
  - Online removal of  $^{222}\text{Rn}$  from cable conduits with cold synthetic charcoal.
- Ability to influence and control liquid xenon convection in the TPC
  - Constrain/reject problematic  $^{222}\text{Rn}$  daughters such as  $^{214}\text{Pb}$ .
- Insert short-lived sources such as  $^{83\text{m}}\text{Kr}$  and  $^{220}\text{Rn}$ .
- Insert and remove via gettering tritiated methane for ER calibration.
- Gas phase Xe recovery
  - Two all-metal gas compressors
  - 144 custom high-purity compressed gas cylinders
  - Diesel generator backup to mitigate power disruptions
- Integrated and intensive purity monitoring campaign with cold-trap mass spectrometry



# TPC Liquid Flow Design



## GAS PHASE AND ELECTROLUMINESCENCE REGION



- Liquid removal from the top of the TPC via weir spill-over mechanism.
  - Sets precise LXe level (~10 micron control)
  - Floating particulates, fibers, dust are collected and removed from the TPC
- However: Xe skin region instrumented as S1-only gamma/neutron veto
  - LXe drain pipe is forbidden in the skin.
  - Drain placed in insulating vacuum space, between inner and outer cryostats.
- Introduce purified LXe via seven inlet tubes integrated into lower PMT array.
- Careful thermal design allows temperature of inlet fluid to be controlled.
  - Inlet liquid temperature allows LXe convection to be promoted or suppressed.



## Weir Return from TPC


TPC Weir  
Port

1/2" Weir  
Drain Line

"Y" fittings


Cooling HX connects to  
Evaporator with thermally  
conductive copper strap  
10W cooling

3/4" Weir Drain Line  
connects to vacuum  
jacketed hose inside  
Special Tee




- Three weir drain located at 120° intervals around the TPC.
- One liquid liter per minute of LXe drains through these ports.
- Vapor lock mitigation:
  - cooling strap supplies 10 W of cooling power to the triple-junction.
  - a dedicated gas return path (not depicted)

# Weir Drain Assembly – TPC insertion into Inner Cryostat




Inner Cryostat Vessel




Flexible weir drain is threaded through port and sealed to a 1/2" SS drain tube which also seals to the ICV vessel.

# LXe Tower: interface to the external Xe handling system




- Purified LXe is condensed and subcooled in the LXe tower. Located outside the LZ water shield at an elevation below the TPC.
- A vacuum insulated conduit connects the bottom of the inner cryostat to the LXe tower.
- Inside the conduit:
  - Two vacuum-insulated cryogenic lines supplying purified LXe to the bottom of the detector.
    - One feeds the TPC active volume
    - The other supplies and flushes the Xe skin region
  - A third cryogenic line returns the LXe from the weir drain to the LXe in preparation for evaporation and gettinger.
- Proportional cryovalves control the flow rates in the three lines.

# LXe Tower Architecture



Gas Xe supply and return



Liquid Xenon Tower

standpipe

- LXe from the weir drain spills into the reservoir.
- Standpipe construction acts as a flow circuit diode:
  - backflow is not possible
  - dust, fibers, particulates, ect., should become trapped
  - Non-volatile impurities (heavy hydrocarbons) should also become trapped.
- Two-phase heat exchanger mates the LXe and gaseous Xe portions of the circuit.
  - Efficiency is 95%
- LXe Subcooler:
  - removes bubbles and froth from the purified LXe
  - Cools the LXe several degrees below saturation
  - LXe trim heaters make final adjustment to LXe temperature prior to supplying it to the TPC.

LXe to and from cryostat


# SLAC System Test (2015-2019): Purification & circulation system prototyping

## SLAC System Test

- 100 kg LXe TPC
- 50 SLPM circulation rate
- Prototyping of LXe tower
  - reservoir
  - Heat exchangers
  - Weir drain
- Key findings regarding weir drain line construction and assembly.
- Design strategies to prevent vapor lock in weir drain and transfer lines.




## LZ @ SURF






# LZ full scale circulation test @ SURF (2020-21)




- Full height test cryostat
  - Necessary to properly replicate the gravity-driven liquid flow between cryostat and LXe tower.
- Complete system exercise:
  - Xe circulation compressors @ 400 SLPM
  - LXe Subcooler & trim heaters
  - Two-phase heat exchanger
- Confirmed that vapor lock mitigations devised in the SLAC System Test remain effective in the full scale LZ hardware.
- Debugging of liquid nitrogen cryogenics systems; reduction of heat leaks.


# Thermal control by design: $^{222}\text{Rn}$ and its daughters are not fully mixed in LZ





# Thermal control by design: $^{222}\text{Rn}$ and its daughters are not fully mixed in LZ



LZ SR1 data:  $^{218}\text{Po}$



LZ SR1 data:  $^{214}\text{Po}$



Rn222 ( $\mu\text{Bq/kg}$ )	Pb214 ( $\mu\text{Bq/kg}$ )	Po214 ( $\mu\text{Bq/kg}$ )
$4.37 \pm 0.31$ (stat)	$3.26 \pm 0.13$ (stat) $\pm 0.57$ (sys)	$2.56 \pm 0.21$ (stat)

# All-metal Diaphragm Compressors


Two Circulation Gas Xe Compressors



Two Recovery Gas Xe Compressors




144 high purity Xe gas cylinders





- Triple diaphragm construction.
- Copper plating on diaphragm forms all-metal seal with compressor heads.
- Check valves sealed to the heads with copper gaskets.

# SAES/Entegris PS5-MGT50 Megatorr Getter



- SAES conservatively recommends two parallel MGT50 units at Xe gas flow rates  $> 275$  SLPM.
- To limit Rn emanation burden from the getter, we use only one.
- We find no degradation in performance up to 500 SLPM.
  - The purifier is able to establish and maintain a temperature of 400 C.

# Electronegative and Methane Removal Performance




- Electron lifetimes in LZ are near 10 ms; far better than the 1 ms requirement.
- To test the getter's ability to remove methane, a 0.07 gram injection was performed. Cold trap mass spec. measurements confirmed good removal with a time constant less than four days.

# Radon Emanation measurements of SAES purifiers / St707 getter pills

- Getter pills have high surface area and operate at 400 C;  $^{222}\text{Rn}$  is expected to be emanated.
- Literature measurements suggest a factor of eight variation in  $^{222}\text{Rn}$  emanation depending on the unit / lot.
- The LZ getter is responsible for about 5% of the  $^{222}\text{Rn}$  in LZ.
  - Literature results suggest this could have been reduced to 1% with the most favorable getter pill lot.

Device	Model	Temp	$^{222}\text{Rn}$ Activity ( $\mu\text{Bq}$ )	Getter Pill Mass (kg)	Activity ( $\mu\text{Bq}/(\text{kg-of-getter})$ )
Maryland [1]	MT50	400 C	142 + 39 - 35	4.4	32±8
Xenon1T (ID #16) [2]	MT50	400 C	1170 ± 150	4.4	265 ± 38
Xenon1T (ID #17) [2]	MT50	400 C	240 ± 30	4.4	55 ± 7
LZ Getter [3]	Megatorr PS5 –MGT50	400 C	2260 ± 2270	15	150 ± 15



St707 Getter Pills




[1] J.E. Armstrong et al., *AIP Conf. Proc.* 2908, 070002 (2023) <https://doi.org/10.1063/5.0161727>

[2] E. Aprile et al. (XENON), *Eur. Phys. J. C* 81, 337 (2021), arXiv:2009.13981

[3] D.S. Akerib et al. (LZ), *Eur. Phys. J. C* (2021) 80:1044, arXiv:2006.02506

# Radon emanation measurements of the purification system & TPC




- Comprehensive Rn screening campaign.
- Most of the Xe circulation components were screened at remote facilities using portable Rn traps.
- The single largest contributor appears to be the inner cryostat.




# Two cold-trap mass spectrometry systems underground at SURF monitor Kr, Ar, N<sub>2</sub>, He, CH<sub>4</sub>

Mobile Sampling System (MSS)




Stationary Sampling System (SSS)




- 77 K cold trap removes most Xe while allowing noble gases and other impurities to pass through.
- Purities are observed with a downstream RGA.
- Sensitivity to Kr is 7 parts-per-quadrillion (g/g)
- Mobile system can be moved and plumbed in to various Xe system components as needed in the Davis cavern. Cold trap is LN cooled.
- Stationary system use a PTR cold head. It is permanently plumbed to the detector.

# Two cold-trap mass spectrometry systems underground at SURF monitor Kr, Ar, N<sub>2</sub>, He, methane

## 77 K Cold Trap




340 ppq (g/g) Kr signal observed by RGA




Note: ppq = parts-per-quadrillion

# Kr measurements in LZ





- LZ requirement was  $< 300$  ppq Kr (g/g) at the start of SR1.
- Achieved Kr concentration was near 100 ppq Kr (g/g).

# Kr measurements in LZ



- Slow rise in Kr concentration is expected due to detector outgassing.
- Rate of increase is 50 ppq (g/g) per year. This meets the LX requirement.
- Observed increase was predicted based upon outgassing constants measured during design phase.

# Charcoal Gas Chromatography for Kr removal at SLAC (2020-2021)





- Kr concentration reduced from  $\sim 10$  ppb to  $123 \pm 10$  ppq (g/g).
- Purified Xe deliveries to SURF complete in September 2021.

Thanks to our sponsors and 35 participating institutions!



U.S. Department of Energy  
Office of Science



Science and  
Technology  
Facilities Council



ibs Institute for  
Basic Science




Questions?



# Liquid System Elevation Plan

Gas Xe supply/return



Three LXe transfer lines inside conduit